



Autogen Bioclear UK Ltd, Holly Ditch Farm, Mile Elm, Calne, Wilts, UK
 Tel: 01249-819008 Fax: 01249-817266 e-mail: info@autogenbioclear.com
www.autogenbioclear.com



Joey™ Gel Casting Systems

Cast gels with ease!

Available separately and included with each Penguin package set, the patented Joey Gel Casting System uses an arrangement of polyethylene bags, or "Gel Pouches", and a simple clamping system to cast polyacrylamide gels. Glass plates and spacers are placed in the pouch and acrylamide is poured between the plates. A Spacer Placer™ is included with each unit to easily align the side spacers. During polymerization, acrylamide is completely contained within the fitted pouch and will not leak. With the Joey Gel Casting system there is no need for tape or a bottom spacer; the convenient pouch design cuts set-up time to minutes. Up to four gels may be cast at one time and pouches may be heat sealed for later use. Twenty-five gel pouches are included with each system.

Model JGC-4

Gel Size: 10cmW x 10cmL

Item Description	Model No.
Complete System	JGC-4

Accessories

Gel pouches, pk of 25	GP4-25
Gel pouches, pk of 100	GP4-100
Spacer Placer, pk of 3	JG4-PL

Model JGC-2

Gel Size: 16cmW x 14, 16cmL

Item Description	Model No.
Complete System	JGC-2

Accessories

Gel pouches, pk of 25	GP2-25
Gel pouches, pk of 100	GP2-100

Model JGC-3

Gel Size: 20cmW x 20cmL

Item Description	Model No.
Complete System	JGC-3

Accessories

Gel pouches, pk of 25	GP3-25
Gel pouches, pk of 100	GP3-100
Spacer Placer, pk of 3	JG3-PL





Autogen Bioclear UK Ltd, Holly Ditch Farm, Mile Elm, Calne, Wilts, UK
Tel: 01249-819008 Fax: 01249-817266 e-mail: info@autogenbioclear.com
www.autogenbioclear.com



Multiple Gradient Caster

Versatile caster for all types of gels

Owl offers the Multiple Gradient Caster for casting gradient or single percentage gels. This leak resistant unit allows you to pump the acrylamide through a port on the base of the caster. This port can be closed to pour gels from the top. Foam and acrylic displacement plates are provided to hold gel plates snugly in the proper orientation and ensure even pressure across the gel. Up to five gels can be poured simultaneously. The Gradient Caster is available in sizes compatible with all Penguin, Puffin and large format single sided systems.

Multiple Gradient Caster

Item Description	Model No.
Multiple Gradient Caster 10cmW x 10cmL	P7-CST
Multiple Gradient Caster 16cmW x 14cmL	P1-CST
Multiple Gradient Caster 20cmW x 20cmL	P2-CST

